

# Pushrims & pushrim covers

Fitting options

Pushrims in aluminium

Coated aluminium pushrims

Pushrims with knobs

Ergonomic pushrims

Steel pushrims

Silicone covers

# Pushrim fitting options

## Fitting options – Double wall rims

### TYPE A:

Pushrim with splices fitted on rim with rivnuts



- » The pushrim can be exchanged without removing the tyre and the tube

## Fitting options – Single wall rims

### TYPE B:

Pushrim with splices fitted with screw and lock nut



- » necessary to remove tyre and tube to exchange pushrim

### TYPE C:


Pushrim with rivnuts fitted with anchor-plates



- » Distance between pushrim and rim is 17 mm as standard. Other lengths on request



**Examples of ergonomic  
pushrims in our Omobic  
offer – see [MBL.DK](http://MBL.DK)**



## Aluminium pushrims with rivnuts black powdercoated, silver anodized or black anodized




- » Ø19 x 1,25
- » 4 or 6 rivnuts
- » Ø19,6 x 1,25
- » Black powdercoated, silver anodized or black anodized

Wheel Size	SD fitting distance	Fit rim types	Weight g	Fitting	Black powdercoated	Silver anodized	Black anodized
20"	426	SW1	-	4 rivnuts	322001	322002	-
22"	464	SW1	-	4 rivnuts	322201	322202	322242
22"	469	DW5	298	6 rivnuts	-	322284	322285
24"	515	SW1	-	4 rivnuts	322401	322402	322442
24"	520	DW5	330	6 rivnuts	-	322495	32240BX
26"	565	SW1	-	4 rivnuts	322601	322602	322642
20"	426	SW1	-	6 rivnuts	322051	322052	-
22"	464	SW1	-	6 rivnuts	322251	322252	-
24"	515	SW1	-	6 rivnuts	322451	322452	322450
25"	534	SW1	-	6 rivnuts	-	322561	322562
25"	539	DW5	341	6 rivnuts	-	322595	322594
26"	565	SW1	-	6 rivnuts	322651	322652	-

# Aluminium pushrims with splices

silver anodized or black anodized



- » Ø19 x 1,25
- » Splices (392115) with two holes
- » Extra long splices used for pushrims which fit extra wide wheels
- » Splices can be adjusted +/- 2 mm

Wheel Size	SD fitting distance mm	Fit rim types	Weight g	Fitting	Splice Type	Silver anodized	Black anodized
20"	422	DW2	-	6 splices	392115	on request	on request
20"	434	SW1	-	4 splices	-	342001	342002
20"	434	SW1	292	6 splices	392115	342051	342052
22"	460	DW2	-	6 splices	392115	on request	on request
22"	466	DW4 & DW7	311	6 splices	392115	342255	342256
22"	472	SW1	315	6 splices	392115	342251	342252
22"	472	SW1	-	4 splices	-	342201	342202
22"	472	SW1	-	6 splices	-	342251	342257
24"	480	Buddy SW Wide	-	6 splices	392144	-	3424E4
24"	510	DW2	338	6 splices	392115	344458	342458
24"	517	DW4 & DW7	343	6 splices	392115	342455	342456
24"	523	SW1	-	4 splices	-	342401	342402
24"	523	SW1	350	6 splices	392115	342451	342452
25"	530	DW2	351	6 splices	392115	342568	342567
25"	542	SW1	359	6 splices	392115	342561	342558
26"	560	DW2	-	6 splices	392115	on request	on request
26"	567	DW4 & DW7	374	6 splices	392115	342655	342656
26"	573	SW1	378	6 splices	392115	342651	342652

# Aluminium pushrims coated for extra grip



## Coated pushrim with splices

- » Ø19 x 1,25
- » Thickness of dipping app 1 mm
- » Fitting with splices (392115) with two holes
- » Splices can be adjusted +/- 2 mm

Wheel size	SD fitting distance mm	Fit rim types	Weight g	Fitting	Standard Grip
20"	428	DW4 & DW7	-	6 splices	342026
20"	434	SW1	-	(check)	342004
22"	466	DW4 & DW7	-	6 splices	342214
22"	472	SW1	-	6 splices	342254
24"	510	DW2	-	6 splices	34249F
24"	517	DW4 & DW7	-	6 splices	342489
24"	523	SW1	-	6 splices	342454
25"	530	DW2	-	6 splices	3425D5



## Coated pushrim with rivnuts

- » Ø19 x 1,25
- » Thickness of dipping app 1 mm
- » Fitting with rivnuts

Wheel size	SD fitting distance mm	Fit rim types	Weight g	Fitting	Standard Grip
20"/451	426	SW1		4	322010
22"/489	464	SW1		4	322205
24"/540	515	SW1		4	322410
26"/590	565	SW1		4	-
20"/451	426	SW1		6	322053
22"/489	464	SW1		6	322260
24"/540	515	SW1		6	322460
25"/559	534	SW1		6	3225D1
26"/590	565	SW1		6	322660

## Aluminium pushrims with knobs




### Aluminium pushrim with 8 knobs and splices

- » Ø24 x 1,25
- » 4 or 6 splices
- » Note: the knobs can be adjusted to different angles
- » Length of knobs: 20 mm
- » Black powder coated

Wheel size	SD fitting distance mm	Fit MBL rim types	Weight (g)	Fitting	Part. no.
20"	435	SW1	-	4 splices	342024
22"	473	SW1	-	6 splices	342274
24"	524	SW1	-	6 splices	342474
26"	576	SW1	-	6 splices	342674

## Ergonomic pushrims



### Ergonomic plastic pushrim with rivnuts

- » Good grip and comfortable to touch also in cold weather
- » Generates heat when braking at high speed
- » Ø19 x 26
- » Non allergic and insulates against cold
- » 4 or 6 rivnuts

Wheel size	SD fitting distance mm	Fits MBL rim type	Weight (g)	Fitting	Part no.
20"	426	SW1	-	4 rivnuts	312002
22"	464	SW1	-	4 rivnuts	312202
24"	515	SW1	-	4 rivnuts	312402
20"	426	SW1	-	6 rivnuts	312052
22"	464	SW1	-	6 rivnuts	312252
24"	515	SW1	-	6 rivnuts	312452
26"	565	SW1	-	6 rivnuts	312652

# Stainless steel pushrims

## brushed or polished



### Stainless steel pushrims with splices

- » Very durable and scratch resistant
- » 6 splices (392169), 2 mm thick
- » Splices can be adjusted +/- 2 mm
- » Ø19 mm

Wheel size	SD fitting distance mm	Fit Omobic rim types	Weight g	Fitting	Polished	Brushed
20"	434	SW1	-	4 splices	352022	352012
22"	466	DW4 & DW7	668	6 splices	352258	on request
22"	472	SW1	-	4 splices	on request	352212
22"	472	SW1	676	6 splices	352272	352262
<b>24"</b>	<b>510</b>	<b>DW2</b>	<b>727</b>	<b>6 splices</b>	<b>352478</b>	on request
24"	517	DW4 & DW7	736	6 splices	352458	on request
24"	523	SW1	-	4 splices	352421	352415
24"	523	SW1	745	6 splices	352472	352462
25"	535	DW4 & DW7	763	6 splices	352558	on request
26"	573	SW1	816	6 splices	352672	352662
26"	567	DW4 & DW7	803	6 splices	352658	on request



### Stainless steel pushrims with rivnuts

- » Very durable and scratch resistant
- » 4 or 6 rivnuts
- » Ø19 mm

Wheel size	SD mm	Fits MBL rim type	Fitting	Polished	Brushed
20"	426	SW1	4 rivnuts	332004	332005
20"	426	SW1	6 rivnuts	332054	332055
22"	464	SW1	6 rivnuts	332254	332205
24"	515	SW1	6 rivnuts	332454	332455
26"	565	SW1	6 rivnuts	332654	332655



## Silicone covers with ribbed or smooth surface



- » Pushrim covers in silicone are a perfect solution for many wheelchair users
- » Can be used onto your existing pushrim to improve grip and comfort
- » Easy to mount in seconds without using any tools
- » These covers are ideal for chipped and damaged pushrims
- » Insulate from cold steel/aluminium
- » Provides better control over a chair as well as better efficiency of propulsion and braking
- » Do not cause allergic reaction
- » Other dimensions are available on request

Pushrim dimension Ød x Ød (mm)	Wheel size	Weight (g)	Part no. black ribbed (1 pcs)	Part no. black smooth (1 pcs)
Ø16-19 x Ø430-437	20"	400	382031	382021
Ø16-19 x Ø438-445	20"	400	382032	382022
Ø16-19 x Ø446-453	20"	400	382033	382023
Ø16-19 x Ø454-461	20"	400	382034	382024
Ø16-19 x Ø462-469	22"	450	382231	382221
Ø16-19 x Ø470-477	22"	450	382232	382222
Ø16-19 x Ø478-485	22"	450	382233	382223
Ø16-19 x Ø486-493	22"	450	382234	382224
Ø16-19 x Ø494-501	22"	450	382235	382225
Ø16-19 x Ø502-510	22"	450	382236	382226
Ø16-19 x Ø511-518	24"	500	382431	382421
Ø16-19 x Ø519-525	24"	500	382432	382422
Ø16-19 x Ø526-532	24"	500	382433	382423
Ø16-19 x Ø533-539	24"	500	382434	382424
Ø16-19 x Ø561-570	26"	550	382631	382621
Ø16-19 x Ø571-580	26"	550	382632	382622
Ø16-19 x Ø581-590	26"	550	382633	382623
MBL - oval plastic	20"	400	-	382030
MBL - oval plastic	22"	450	-	382230
MBL - oval plastic	24"	500	-	382430