



Trays




Tray in transparent plastic – with or without mounting system

- » Table material: PETG
- » High surface hardness
- » Impact resistant
- » Shatter-proof
- » Mounting system: Black KTL steel
- » Easy adjustment of width



Description	Part no.
Complete tray incl. mounting system for 320–410 mm seat width	721059
Complete tray incl. mounting system for 460–510 mm seat width	72105J
Tray excl. mounting system for 320–410 mm seat width	721004-1
Tray excl. mounting system for 460–510 mm seat width	721006-1



Mounting system for trays

- » Black KTL steel
- » Easy adjustment of width

Description	Part no.
Complete mounting system for 320–310 mm seat width	72114LR
Complete mounting system for 460–510 mm seat width	72114LR