

Aluminium pushrims Coated for extra grip

Standard grip (893201)

- Black foam 65 SHA
- Based on PVC
- Thickness 1 mm
- Fulfill all EU regulations
- 1,1 friction coefficient

Lizard grip – Extra sticky (893204)

- Extra sticky – good for sport
- Black foam 65 SHA
- Based on PVC
- Thickness 1 mm
- Fulfill all EU regulations
- 0,7 friction coefficient



- Ø19 x 1,25
- Thickness of dipping app 1 mm
- Fitting with 6 splices or rivnuts
- Splices (392115) with two holes
- Splices can be adjusted +/- 2 mm

Wheel Size	SD / rivnuts fitting distance (mm)	Fit rim types	Weight g	Fitting	Standard Grip	Lizard Grip Produced in Poland
20"	428	DW4 & DW7	480	6 splices	342026	-
22"	466	DW4 & DW7	520	6 splices	342214	-
24"	510	DW2	565	6 splices	34249F	-
24"	517	DW4 & DW7	570	6 splices	342489	3424A5
24"	523	SW1	580	6 splices	342454	3424H7
25"	530	DW2	585	6 splices	3425D5	-
24"	515	SW1	570	6 rivnuts	322460	-

Other sizes on request

